


Post-Operative Care



The kit includes a separate envelope labeled "NOT FOR USE IN O.R." that includes a **standard 3 cc Luer Lock irrigating syringe and adapter tip**. During the post-op period, irrigating the airway tubes will relieve any clogging of mucus or blood. The enclosed **PATIENT HOME CARE INSTRUCTIONS** card should be reviewed by the recovery room nurse with the patient and caregiver.

Office Removal Instructions

Prior to removing the airway, septal splints, and packing if present, instill a topical anesthetic/decongestant mixture into the nasal passages. Shrinking and anesthetizing the mucosa will facilitate a rapid, pain-free removal and thus satisfactory patient experience. Using a hemostat clamp, grasp the wall of each tube and extract the device.

If the Reltok Ultra-Smooth Septal Splint has been employed, remove the airway first for easier access to the splint. Further topically anesthetize and decongest the septum prior to splint release and removal.



7740 Records Street • Indianapolis, IN
46226 800-428-1610 • 317-545-6196
www.anthonypproducts.com

anthonypproducts.com/reltok


The Reltok Clear-Flo™ Nasal Airway is FDA cleared.
Clear-Flo Nasal Airway, and Always Clear, Always Comfortable are trademarks of Reltok Nasal Products & Anthony Products, Inc.

©2025 Reltok Nasal Products, LLC
U.S. Patent No. 8,092,478 and 8,874,486

RELTOK CLEAR-FLO™
NASAL AIRWAY

The Reltok Clear-Flo™ Nasal Airway

ALWAYS CLEAR, ALWAYS COMFORTABLE™



Instruction For Use

The Reltok Clear-Flo™ Nasal Airway System kit consists of:

- The Reltok Clear-Flo™ Nasal Airway, made of medical grade silicone, latex-free.
- The Reltok Ultra-Smooth Septal Splints, Class VI medical grade silicone, Parylene coated enabling a friction-free insertion.
- A standard flexible 10Fr suction catheter.
- A 3 cc Luer Lock irrigating syringe and blue adapter tip.

The septal splints, nasal airway and suction catheter are for use during surgery. The irrigating syringe and tip are for post-operative irrigation.

INDICATIONS FOR USE:

The Reltok Ultra-Smooth Septal Splints are designed to stabilize the reconstructed nasal septum and prevent adhesions between the septum and other intranasal structures.

Reltok Clear-Flo™ Nasal Airway is used to provide nasal/sinus surgery patients' immediate post-operative period with satisfactory nasal air flow, whether or not the nose is packed or has had splints or stent placed within the nasal cavity. The device provides the anesthesia specialist unhindered access for pharyngeal and supra-laryngeal suction following completion of the surgical procedure, in the operating room or recovery room.

The Reltok Clear-Flo™ Nasal Airway is used to provide nasal/sinus surgery patients' immediate post-operative period with satisfactory nasal air flow, whether or not the nose is packed or has had splints or stent placed within the nasal cavity. The device provides the anesthesia specialist unhindered access for pharyngeal and supra-laryngeal suction following completion of the surgical procedure, in the operating room or recovery room.

WARNING: As with any surgical procedure, care should be exercised in the insertion and maintenance of these devices.


PRECAUTIONS: Inspect package before use. DO NOT USE if package has been opened or is damaged. The contents of this package are for single use only. DO NOT re-sterilize. DO NOT use after expiration date.

CAUTION: Federal law restricts this device for sale by, or on the order of, a license physician.

Surgeon Insertion Instructions

The Reltok Clear-Flo™ Nasal Airway is introduced at the conclusion of the operation before inserting any nasal packing or gels. For ease of insertion, lubricate the airway tubes with saline, ointment, or lubricating jelly. After initially inserting the airway partly into the nasal cavity, use a standard, thin-tip nasal speculum to inspect the nasal cavity to ascertain the position of the tubes.

Using direct vision, advance the airway further into the nose. Next, use the inferior portion of the speculum blade or bayonet forceps to direct each airway tube downward onto the nasal floor.





The airway tube will snap into place on the floor of the nose and maintain that position lateral to the pre-maxillary bone and medial to the inferior turbinate.

When both airway tubes are properly seated, the bridge connecting them will be flush against the columella. **NOTE: If an open rhinoplasty procedure has been performed,** the surgeon may wish to divide the bridge and secure each airway tube separately, rather than have the bridge contact the transcolumellar incision.



After inserting and seating the airway, pass the **10 Fr plastic suction catheter** through each tube and suction fluid from the pharynx. This maneuver also confirms that the back opening of the device is unobstructed. Later, the anesthesia specialist can use the same flexible suction catheter to remove blood and mucus from the throat.



While the packing of choice or gel is placed, the nasal speculum stabilizes the airway tube.